

The diagram is structured as follows:

- Top Level:** An umbrella shape with a black outline and a blue oval in the center containing the text "SERVICE LIFE PREDICTION OF POLYMERIC SYSTEMS".
- Second Level:** A light blue rectangular box with a black border containing the text "Contributions from Individual Constituents" in green.
- Third Level:** A grey upward-pointing arrow connects the second level to a red horizontal line.
- Third Level:** A red horizontal line with the text "Metrology/Modeling" in red, underlined, centered above it.
- Bottom Level:** Three light blue rectangular boxes with black borders, each connected to the red line by a red arrow pointing downwards.
  - Left Box:** Contains the text "Photoreactivity of TiO<sub>2</sub> Nanostructures".
  - Middle Box:** Contains the text "Multi-scale Structure and Dispersion Measurements".
  - Right Box:** Contains the text "Surface Mechanical Characterization". This box is highlighted with a thick red border.

**SERVICE LIFE PREDICTION  
OF POLYMERIC SYSTEMS**

**Contributions from  
Individual Constituents**

**Metrology/Modeling**

Photoreactivity  
of TiO<sub>2</sub>  
Nanostructures

Multi-scale Structure  
and Dispersion  
Measurements

Surface  
Mechanical  
Characterization





# Mechanical Property Characterization of Polymer Film Surfaces and Interfaces


---

Aaron M. Forster




# Constituent Contribution to Service Life Prediction

New Technical Idea





# What is the Problem?



- Dispersion, Photoreactivity, Matrix Chemistry (**micro and nano-scale**) contribute to the chemical, mechanical, and physical changes that are associated with macroscale coating degradation



# Why Nanoindentation?




---

- Difficult to accurately measure the mechanical properties of small volumes at the surface of polymer films
  - **Nanoindentation**
  - Brillouin Light Scattering
  - Quartz Crystal Microbalance
  - Other Acoustic Methods
- **Nanoindentation or Instrumented indentation testing (IIT)** can provide mechanical property data at length scales that are several orders of magnitude less than bulk. micron size spatial resolution across exposed surface
  - measure 1  $\mu\text{m}$  to 3  $\mu\text{m}$  into surface
  - quantify a range of elastic and viscoelastic properties





# Principles of Indentation



MTS Nano  
XP/CSM


Travel = 2 mm  
 Displ. Res. = < 0.01 nm  
 Max Load = 1 kg  
 Load Res. = 50 nN





# Quasi-static / Dynamic

## Continuous Stiffness Method



$h_{\max}$  = maximum displacement  
 $h_f$  = final depth  
 $h_i$  = intercept displacement

$$E_r = \frac{\sqrt{\pi}}{2} \left( \frac{S}{\sqrt{A}} \right)$$

- Dynamic oscillation is superposed over a given loading history
- Settings:
  - harmonic amplitude = 1-50 nm
  - frequency target = 10-250 Hz
- Better sensitivity to surface contact, provides continuous estimate of E

## Dynamic Mechanical Characterization

Oscillate tip at sample surface

$$E_r' = \left[ \frac{P_0}{\Delta h_0} \cos \delta \right] \frac{\sqrt{\pi}}{2\sqrt{A}} = \frac{S \sqrt{\pi}}{2\sqrt{A}}, \text{ elastic}$$


$$E_r'' = \left[ \frac{P_0}{\Delta h_0} \sin \delta \right] \frac{\sqrt{\pi}}{2\sqrt{A}} = \frac{C \omega \sqrt{\pi}}{2\sqrt{A}}, \text{ viscous}$$





# Technical Approach

## *Incorporation of Surface Mechanical Property Measurements*






# Why is this difficult?

---

## *Incorporation of Surface Mechanical Property Measurements*



- Contact radii are 2 – 10 micron and depths of 1 – 3 micron
- Tip-sample interaction must be characterized
  - Roughness
  - Presence of particles and particle geometry
  - Effects of tip geometry
- Multiple approaches for polymer analysis
  - No stand-outs
  - No clear method for difficult samples

Relationship between dispersion, chemistry, and  
mechanical properties for small volumes





# Gantt chart - Organization

Y1

Y2

Y3

Y4

## Quasi-static measurements

Modulus and Hardness – comparison to bulk

Link to failure indicators

Link to particle, dispersion, surface characteristics

## Dynamic Measurements

Creep and Stress Relaxation Dynamic Mechanical (E' and E'')

Comparison to bulk

Link to failure indicators and film characteristics

Continually refine and improve current measurement techniques and advance data analysis

## New Techniques\*



New metrologies for increasing resolution and accuracy



\*: investigate additional mechanical property measurement methods



# Particle Dispersion Effects Modulus



- Quantifying dispersion is important
- Solid line is the theoretical model  $\longrightarrow E_r = 1 + 2.5\phi + 14.1\phi^2$
- Particle volume concentration ( $\phi$ ) is not accurate with nanoparticles


\*Error bars represent one standard deviation






# Particle Dispersion in Films - SEM


Acrylic Urethane  
Filler B (250 nm), 2.5%



Acrylic Urethane  
Filler C (20 nm), 2.5%




Epoxy  
Filler B (250 nm), 2.5%



- ☐ SEM image of sectioned films (microtomed) before degradation
- ☐ The lack of dispersion in AU Filler C and Epoxy Filler B are evident
- ☐ AU Filler B has much better dispersion, but clumps are still present



# No Difference in Pigment for Filled Epoxies




\*Error bars represent 95% confidence levels





# Differences in Pigment for Acrylic Urethane





*Are the differences dispersion or photoreactivity?*

\*Error bars represent one standard deviation





# Styrene-Acrylic Latex Exhibits Less Stiffening



<sup>†</sup>AFM courtesy of X. Gu

*Difference in Film Morphology for Latex*


\*Error bars represent one standard deviation






# Elastic Modulus Quantified Across the Interface

## Nanoindentation



## Confocal



X. Gu, P. Drzal, et al.

Fresh Epoxy: Modulus (3.8 Gpa  $\pm$  0.1Gpa)



Degraded Sample: Surface Modulus (5.0 Gpa  $\pm$  0.1Gpa)





# Dynamic Measurements depend on Film Properties

Glassy, Stiff polymers



- Tip-sample interactions
- Viscoelastic materials


P. Drzal, M. VanLandingham, et al.





# Storage Modulus Increases with Crosslinking

Bulk mechanical measurements of a rubbery sealant material exposed to accelerated weathering conditions



❑ Dynamic mechanical thermal analysis: Temperature sweep undertaken at 1Hz from -100 $^{\circ}\text{C}$  to 120 $^{\circ}\text{C}$ .

❑ Indication for **macroscopic cross-linking** occurs.

❑ Transition between glassy modulus and rubbery modulus become less distinctive, indicating specimens become more **heterogeneous**.

Exposure environment: 30 $^{\circ}\text{C}$  and 0% relative humidity

\*C. White, NIST





# Summary

---

- Conduct modulus measurements on filled/unfilled coatings
  - Epoxy, acrylic urethane, and latex
  - Bulk measurements on sealant materials
- Track changes in modulus at the surface of the coatings throughout the degradation cycle
  - Epoxy and acrylic urethane
- Challenges remain:
  - Effect of roughness
  - Influence of particle size, dispersion
  - Viscoelastic characterization of degraded surfaces
  - Tip-Sample interaction





# Gantt chart – Progress Metrics

---

Y1

Y2

Y3

Y4

## Quasi-static measurements

---

- ❑ Modulus – comparison to bulk
- ❑ failure indicators – identified
- ❑ correlations to dispersion, photoreactivity, and chemistry

## Viscoelastic Measurements

---

- ❑ Storage and loss modulus – ability to measure and comparison to bulk
- ❑ Creep and stress relaxation - ability to measure and comparison to bulk
- ❑ Correlations to dispersion, photoreactivity, and chemistry
- ❑ Modeling the tip-sample interaction





# Impact

---

- Develop methodologies that **relate** mechanical properties (quasi-static and dynamic) to durability analysis
- A **link** between filler chemistry, dispersion, and mechanical properties that is incorporated into service life prediction
  - Scales from constituent to coating performance
- Fundamental material system chosen and evaluated for **performance metrics**
  - engineered durability – cheaper, better, faster!





# Acknowledgements

---

- ❑ Pete Drzal – PPG
- ❑ Stephanie Watson – NIST
- ❑ Amanda Forster - NIST
- ❑ LiPiin Sung – NIST
- ❑ Xiaohong Gu – NIST
- ❑ Jeff Comer – Cal Poly
- ❑ Laura Johnson – Cal Poly
- ❑ Tom Juliano – Army Research Laboratory







# The End


---

Questions?




# Particle Dispersion in Films - Confocal


## Acrylic Urethane - 2.5 % Filler C



Pigment Clumps: poor dispersion




Holes!






# Film Topography After Degradation

## Confocal Microscopy



$Sq = 0.067 \mu\text{m}$





$Sq = 0.181 \mu\text{m}$

- ☐ 2D projection image (150x), scale bar is 20 μm, area is 60 μm x 60 μm
- ☐ Increase in surface roughness, exposure of pigment, and pitting occur with increased exposure (20 weeks)
- ☐ Presents a challenge for indentation and ATR measurements




# Dynamic Mechanical Characterization



$$E'_r = \left[ \frac{P_0}{\Delta h_0} \cos \delta \right] \frac{\sqrt{\pi}}{2\sqrt{A}} = \frac{S \sqrt{\pi}}{2\sqrt{A}}$$

$$E''_r = \left[ \frac{P_0}{\Delta h_0} \sin \delta \right] \frac{\sqrt{\pi}}{2\sqrt{A}} = \frac{C \omega \sqrt{\pi}}{2\sqrt{A}}$$



Note:


- Works well on glassy and elastomeric materials
- Sample volume is changing
- Assumes harmonic unloading is elastic
- Most polymers are viscoelastic!
- Soft materials require large contacts
- Beware of treating like a black box






# Dynamic Measurements

Epoxy



PMMA



The agreement is pretty good  
The Epoxy has more strain dependence.





# Dynamic Properties

---

$$E'_r = \frac{S\sqrt{\pi}}{2\sqrt{A}}$$

$$E' = \frac{\sqrt{\pi}}{2\sqrt{A}} \left( \frac{P_0}{h_0} \cos \delta \right)$$


$$E''_r = \frac{\omega C_s \sqrt{\pi}}{2\sqrt{A}}$$

$$E'' = \frac{\sqrt{\pi}}{2\sqrt{A}} \left( \frac{P_0}{h_0} \sin \delta \right)$$

This analysis is a good starting place, but should fail when the materials show any viscoelastic nature. The geometry of the nanoindenter does not approximate the gap loading conditions.









# Motivation

---


- ❑ Commercial coatings add particles (pigments) as filler
  - Improves optical properties
    - ❑ Whiteness, brightness, opacity, ...
  - Improves mechanical properties
    - ❑ Scratch resistance, modulus, ...
- ❑  $\text{TiO}_2$  is a semiconducting material that is commonly used
  - Photoreactivity implications for coating performance
    - ❑ Changes in photodegradation mechanisms
- ❑ Questions:
  - *Does photoreactivity affect weatherability?*
  - *How do surface mechanical properties change with weathering?*
  - *Do mechanical changes correlate to thermal or chemical measurements at the degraded film surface?*






# New Technical Idea

---



Indoor



## Service Life Prediction Clear Coatings

- *indoor now drives  
performance metrics because  
moved to a dosage model*



Outdoor

